



SOLVENTS / PENETRANTS

Gram Decolouriser

DESCRIPTION

An ethanol and acetone blend for laboratory cell staining.

FEATURES

- Fast acting.
- Helps dissolve the lipid layer from the gram negative cells thereby allowing the leaching of the stain into the cells.
- Cost-effective.
- Ready-to-use.

PACKAGING

Gram Decolouriser is available in drums.

SAFETY

Please consult the MSDS on **Gram Decolouriser** before use.

The information and recommendations contained herein are based on data believed to be correct at this time. However, no guarantee of any kind, expressed or implied is made with respect to the information contained herein.

Typical Characteristics

Appearance.....	colourless clear liquid
Odour.....	characteristic
pH.....	n/a
Vapour Pressure (mmHg @ 20°C)	186
Vapour Density (air = 1).....	2 @ 15°C
Boiling Point (°C).....	56-78
Freezing / Melting Point (°C).....	< -100
Flash Point (closed cup)	-18°C

 MSDS available

QUEENSLAND
Tel +61 7 3308 5200
Fax +61 7 3308 5201

NEW SOUTH WALES
Tel +61 2 8709 9600
Fax +61 2 9791 5066

VICTORIA
Tel +61 3 9408 4144
Fax +61 3 9408 4155

SOUTH AUSTRALIA
Tel +61 8 8222 1333
Fax +61 8 8260 7883

WESTERN AUSTRALIA
Tel +61 8 9352 3000
Fax +61 8 9352 3033